

PREempt[®]

Ready to Use One-Step
Surface Cleaner and
Disinfectant formulated

with **ACCELERATED[®]**
HYDROGEN PEROXIDE
Patented



Ready to Use

WHY CHOOSE PREempt™?

- **CLEANER** Effective disinfectant as well as a cleaner resulting in added confidence that disinfection can occur.
- **FASTER** Rapid and realistic contact times on a broad spectrum of difficult to kill pathogens.
- **RESPONSIBLE** Designed to be easier on surfaces, staff and occupants. Users are not exposed to VOCs (volatile organic compounds) and PPE** is not required as per health and safety rating.
- **SUSTAINABLE** The active ingredient, Hydrogen Peroxide, breaks down into water and oxygen leaving no active residues which helps to reduce environmental impact.

DIRECTIONS FOR USE

Disinfection of non-porous hard surfaces such as the exterior of fume hoods, tabletops and other environmental surfaces.

1. Heavily soiled surfaces require cleaning prior to disinfection.
2. Apply to surface with cloth or disposable wipe.
3. *Ensure surface remains wet for 3 minutes.*

Disinfection of Non-Critical, Pre-cleaned Instruments and Medical Devices:

1. Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and dried.
2. Clean and rinse lumens of hollow instruments before filling with this product.
3. Using either tray or an ultrasonic unit, immerse instrument in PREempt™ RTU for 3 minutes.
4. Wipe with a clean, damp cloth or paper towel and allow to air dry.
5. Discard solution after each use.

For use in laboratories, clean rooms and other critical environments that require cleaning and surface disinfection. This includes work stations, fume hoods, laboratory counter tops, equipment and other hard non-porous environmental surfaces.



AVAILABLE PRODUCTS

PREempt™ Ready to Use One-Step Surface Cleaner and Disinfectant

DIN: 02436698

SHELF LIFE: 2 Years

CONCENTRATION OF ACTIVE:
0.5% Hydrogen Peroxide

- **Format Size:** 1 L Bottle
Product Code: 11101
Packaging Specifications: 12x1 L
- **Format Size:** 4 L Bottle
Product Code: 11105
Packaging Specifications: 4x4 L

EFFECTIVE AGAINST A BROAD-SPECTRUM OF BACTERIA, VIRUSES AND FUNGI IN 3 MINUTES!

Virucidal: 3 Min.

- Proven effectiveness against the Poliovirus Type 1, Sabin strain type 1 which allows for a Broad-Spectrum Virucide claim against most enveloped and non-enveloped viruses
- HIV-1 Human Immunodeficiency Virus (HIV)
- Human Coronavirus 229E
- Feline Calicivirus, F9 Strain, as a surrogate for Norwalk and Norwalk-like viruses
- This product has demonstrated effective against Poliovirus and is expected to inactivate all Influenza A viruses including 2009 (H1N1) pandemic Influenza A virus

Bactericidal: 3 Min.

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- Methicillin-resistant *Staphylococcus Aureus* (MRSA)
- Vancomycin-resistant *Enterococcus faecalis* (VRE)

Fungicidal: 3 Min.

- *Trichophyton mentagrophytes*

Tuberculocidal: 3 Min.

- *Mycobacterium terrae*

Broad-Spectrum Sanitizing: 30 Sec.

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Klebsiella pneumoniae*
- *Escherichia coli* O157:H7
- *Campylobacter jejuni*



PREempt™ Ready to Use One-Step Surface Cleaner and Disinfectant is designed for use as a daily cleaner and disinfectant which has easy to use formats with no mixing or measuring required.

Please refer to reference sheet for a complete list of pathogenic organisms, additional features, and claims.

PREempt™ branded products are manufactured by Virox Technologies Inc.

Accelerated Hydrogen Peroxide® (AHP®) and PREempt™ are registered trademarks of Virox Technologies Inc.



For more information on this or our other products, visit:

> virox.com

2770 Coventry Road, Oakville, ON L6H 6R1 CANADA
1-800-387-7578 | www.virox.com



Engineering Revolutionary Disinfectants for the War Against Microbes